

## Connect Your CareStream™ to Lumident's Workflow



Save time and money when you send scans directly from your CareStream™ digital impression scanner. **In addition to all-ceramic restorations, Lumident accepts digital impressions for PFM and full-cast restorations.**

- Use the email [digital@lumident.net](mailto:digital@lumident.net) for your first case.
- Once the initial case is sent, the doctor to lab communication is complete.
- In the comment box, be sure to include the patient name, material to be used and due date.

If you have any questions, contact Lumident at 866-586-4336 and ask for digital impression technical support.

### Looking for more information on our products?

For material descriptions, preparation and cementation guides and billing codes, visit [lumident.net/products](http://lumident.net/products)

Save money on every digital submission by choosing model free when possible.

## SAVE \$20 On Your First Digital Case

Digital submissions are simple when you choose Lumident. Reference coupon code **GO\_DIGITAL** on your Rx when you upload your first case and save \$20 on your first digital submission.

**Coupon Code: GO\_DIGITAL**